



To whom it may concern

Im Wiesengrund 13
D 46342 Velen
Tel.: 02863/95024
info@schaefers-medizintechnik.de
www.schaefers-medizintechnik.de

Geschäftsführerin:
Christa Schaefers

St.-Nr.: 307/5908/2066
Ust.-IdNr.: DE124218493
Amtsgericht Coesfeld; HRB 21828

03.02.2025

Technical description of the medical device

irrimed classic is an electric irrigation pump as a medical device for bowel irrigation or transanal irrigation (TAI). It has a battery-operated pump with which the water pressure can be continuously adjusted electronically. The LI-Ion battery is rechargeable and maintenance-free. It is hermetically sealed and can be installed and operated in any position without electrolyte leakage. It has an energy of 19.24 Wh.

The charger for the battery is a plug-in AC/DC adapter, input 220 V - 240 V , 50 Hz or 60 Hz, output 8.4 V 1.3 A; DC .

The water temperature is shown on the LCD display.

This medical device conforms to Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

The EN standards and statutory regulations are complied with, e.g:

EN 980, EN 1041, EN ISO 14971, EN ISO 10993, EN ISO 60601-1-2, 1-8, 1-11.

The company is certified according to EN ISO 13485:2016

for the manufacturer


Christa Schaefers